

85907007

Monopolar Electrodes for MIS with Ceramic



Technical specifications

Dimensions and Weights

|                    |        |
|--------------------|--------|
| Total length       | 360 mm |
| Total length in cm | 36 cm  |
| Diameter           | 5 mm   |

Basic Specifications

Cleaning and Reprocessing

|                         |             |
|-------------------------|-------------|
| Reusable                | Yes         |
| Washing machine         | 95 °C       |
| Autoclave sterilization | max. 137 °C |

Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | N/A       |
| FDA Approval             | N/A       |

Additional Product Information

|                           |                |
|---------------------------|----------------|
| Accessories rated voltage | max. 2 kVp     |
| Packaging unit            | 1              |
| GTIN 14                   | 04250418752411 |