

89505090

Monopolar Electrode Extension



Technical specifications

Dimensions and Weights

Total length	90 mm
Dimensions	mm

Basic Specifications

Shaft diameter	2,4 mm
----------------	--------

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 4 kVp
Packaging unit	1
GTIN 14	04250418729079