


82650003

Monopolar Coagulation Forceps, Electrode



Technical specifications

Dimensions and Weights

Total length	330 mm
Diameter	5 mm

Basic Specifications

Electrode shape	MIXter Grasping Forceps 90°, cross serrated
-----------------	---

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 2 kVp
Packaging unit	1
GTIN 14	04250418746731