

85902100

Monopolar Micro Needle Electrode



Technical specifications

Dimensions and Weights

Diameter	5 mm
----------	------

Basic Specifications

Electrode shape	Needle
-----------------	--------

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 2 kVp
Packaging unit	1
GTIN 14	04250418735902