


89621110

Monopolar Forceps, straight



Technical specifications

Dimensions and Weights

Total length	210 mm
Total length in inches	8,27 "
Tip width	1,6 mm

Basic Specifications

Forceps shape	straight
---------------	----------

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 2 kVp
Packaging unit	1
GTIN 14	04250418729529