


89714112

Monopolar Forceps, straight



## Technical specifications

### Dimensions and Weights

|                        |        |
|------------------------|--------|
| Total length           | 145 mm |
| Total length in inches | 5,71 " |
| Tip width              | 1,4 mm |

### Basic Specifications

|               |          |
|---------------|----------|
| Forceps shape | straight |
|---------------|----------|

### Cleaning and Reprocessing

|                         |                     |
|-------------------------|---------------------|
| Reusable                | Yes                 |
| Manual reprocessing     | Yes, related to IFU |
| Washing machine         | 95 °C               |
| Autoclave sterilization | max. 137 °C         |

### Approvals

|                          |           |
|--------------------------|-----------|
| Directive 93/42 EEC      | Class IIb |
| Regulation (EU) 2017/745 | N/A       |
| FDA Approval             | N/A       |

### Additional Product Information

|                           |                |
|---------------------------|----------------|
| Accessories rated voltage | max. 2 kVp     |
| Packaging unit            | 1              |
| GTIN 14                   | 04250418729963 |