


89500001

Monopolar Tungsten Needle Electrode, 4mm



Technical specifications

Dimensions and Weights

Total length	64 mm
Dimensions	mm
Diameter	0,5 mm

Basic Specifications

Shaft diameter	4 mm
Electrode shape	Needle

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 4 kVp
Packaging unit	1
GTIN 14	04250418728300