

89217040

## Monopolar Electrode Handle



## Technical specifications

### Dimensions and Weights

Dimensions	mm
------------	----

### Basic Specifications

Shaft diameter	4 mm
----------------	------

### Cleaning and Reprocessing

Reusable	Yes
----------	-----

Manual reprocessing	Yes, related to IFU
---------------------	---------------------

Washing machine	95 °C
-----------------	-------

Autoclave sterilization	max. 137 °C
-------------------------	-------------

### Approvals

Directive 93/42 EEC	Class IIb
---------------------	-----------

Regulation (EU) 2017/745	N/A
--------------------------	-----

FDA Approval	N/A
--------------	-----

### Additional Product Information

Accessories rated voltage	max. 5 kVp
---------------------------	------------

Packaging unit	1
----------------	---

GTIN 14	04250418728171
---------	----------------