

89221040

Monopolar Electrode Handle



Technical specifications

Dimensions and Weights

Dimensions	mm
------------	----

Basic Specifications

Shaft diameter	4 mm
----------------	------

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Accessories rated voltage	max. 5 kVp
Packaging unit	1
GTIN 14	04250418739979