


80100242

Monopolar Cable



Technical specifications

Dimensions and Weights

| | |
|--------------|-----|
| Cable length | 3 m |
|--------------|-----|

Basic Specifications

| | |
|---------------------------------|--------------------------|
| Cable connector instrument side | 4 mm Socket angled |
| Cable plug generator side | Bovie, Valleylab, Commed |

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|----------|
| Directive 93/42 EEC | N/A |
| Regulation (EU) 2017/745 | Class I |
| FDA Approval | Class II |

Additional Product Information

| | |
|---------------------------|----------------|
| Accessories rated voltage | max. 5 kVp |
| Packaging unit | 1 |
| GTIN 14 | 04250418701211 |