


89101072

Monopolar Cable for Electrode Pencil



Technical specifications

Dimensions and Weights

Cable length	4 m
--------------	-----

Basic Specifications

Cable connector instrument side	for Monopolare Electrode Handle
Cable plug generator side	Martin, Berchtold, Aesulap GK 55, GK 60

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	N/A
Regulation (EU) 2017/745	Class I
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 5 kVp
Packaging unit	1
GTIN 14	04250418728027