


82400019

POWERGRIP / ORBITARIS Bipolar Coagulation Forceps, Electrode



Technical specifications

Dimensions and Weights

| | |
|--------------|--------|
| Total length | 340 mm |
| Diameter | 5 mm |

Basic Specifications

| | |
|-----------------|----------------------------|
| Electrode shape | Dissector Maryland, curved |
|-----------------|----------------------------|

Cleaning and Reprocessing

| | |
|-------------------------|---------------------|
| Reusable | Yes |
| Manual reprocessing | Yes, related to IFU |
| Washing machine | 95 °C |
| Autoclave sterilization | max. 137 °C |

Approvals

| | |
|--------------------------|-----------|
| Directive 93/42 EEC | Class IIb |
| Regulation (EU) 2017/745 | N/A |
| FDA Approval | Class II |

Additional Product Information

| | |
|---------------------------|----------------|
| Accessories rated voltage | max. 500 Vp |
| Packaging unit | 1 |
| GTIN 14 | 04250418715195 |