


82408015

POWERGRIP / ORBITARIS Bipolar Coagulation Forceps, Electrode



Technical specifications

Dimensions and Weights

Total length	340 mm
Diameter	5 mm

Basic Specifications

Electrode shape	Grasping Forceps, fenestrated, curved
-----------------	---------------------------------------

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418732932