


82460324

POWERGRIP 3.0 micro Bipolar Coagulation Forceps, Electrode



Technical specifications

Dimensions and Weights

Total length	290 mm
Diameter	3 mm

Basic Specifications

Electrode shape	Scissor, curved
-----------------	-----------------

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Accessories rated voltage	max. 300 Vp
Packaging unit	1
GTIN 14	04250418748216