



86710038

Shafts for Resectoscopy, Exterior Continuous Flow Sheath with 2 Stopcocks rotatable



Technical specifications

Dimensions and Weights

Schaftdurchmesser in Charr	26 Charr
----------------------------	----------

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	N/A

Additional Product Information

Packaging unit	1
GTIN 14	04250418741354