


85530080

## SLIMLINE Bipolar Coagulation Forceps, Electrode



### Technical specifications

#### Dimensions and Weights

Total length	250 mm
Diameter	3,4 mm

#### Basic Specifications

Electrode shape	Rippled Bar
-----------------	-------------

#### Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

#### Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

#### Additional Product Information

Accessories rated voltage	max. 500 Vp
Packaging unit	1
GTIN 14	04250418723824