


85530005

SLIMLINE Bipolar Coagulation Forceps, Shaft



Technical specifications

Dimensions and Weights

Total length	200 mm
Diameter	3,4 mm

Cleaning and Reprocessing

Reusable	Yes
Manual reprocessing	Yes, related to IFU
Washing machine	95 °C
Autoclave sterilization	max. 137 °C

Approvals

Directive 93/42 EEC	Class IIb
Regulation (EU) 2017/745	N/A
FDA Approval	Class II

Additional Product Information

Packaging unit	1
GTIN 14	04250418723695